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- Before handling the Sample Collection Card, thoroughly clean the work surface. With gloved hands, carefully open the Sample Collection Kit and remove the Sample Collection Card. **IMPORTANT:** Wear gloves when handling the Sample Collection Cards to avoid cross-contamination with your own DNA.
- Clearly label the Sample Collection Card with the cell line name/designation. If sending multiple cell lines use a separate card for each cell line and make sure the information on the card matches the information on the Sample Submission Form.
- Spot 40µl of the cell suspension prepared in step 2 above at  $1 \times 10^6$  cells/ml in the center of the circle on the inside of the Sample Collection Card.
- Allow the Sample Collection Card to air dry in a laminar flow hood at room temperature.
- When the Sample Collection Card is dry, place it and one desiccant pack (provided with the Kit) in the Multibarrier Pouch with the barcode visible on the see-through side of the pouch. **\*\*To avoid cross-contamination, use one Multibarrier Pouch per sample being submitted.**
- Be sure to completely close the Multibarrier Pouch to preserve and protect the sample.
- Repeat this process for each sample being submitted for testing, manipulating only one cell line at a time to avoid cross-contamination.
- When the cell lines have been spotted and sealed in separate Multibarrier Pouches, place them and the completed corresponding Sample Submission Forms into the pre-addressed Return Envelope(s). If space allows, you may place multiple Multibarrier Pouches into one Return Envelope (be sure to include the Sample Submission Form(s)).
- Affix appropriate postage and place sealed, pre-addressed Return Envelope(s) in the mail. If using an overnight service, please send to the address listed on the pre-addressed Return Envelope. **IMPORTANT:** For International Customers, please add the words "FTA Sample Papers" to the item description on your Air bill labels to avoid delays when mailing your samples.

- Checklist**  
 Before mailing your sample, did you:
- Complete the Submission Form
  - Spot the cells in the Sample Card and place the card in the Ultra-barrier Pouch
  - Include the Submission Form and Ultra-barrier Pouch containing FTA Sample Card in the pre-addressed Return Envelope
  - If mailing more than one sample, please be sure to match the FTA Card barcode number with the Submission Form barcode number
  - Mail the pre-addressed Return Envelope
  - Add words "FTA Sample Papers" to item description on Airway bill labels
- (Continue on back) →

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**Cell Line Authentication Service**  
STR Profile Report

Sample Submitted By: [Redacted]

Email Address: [Redacted]

ATCC Sales Order: SOJ25378  
FTA Barcode: STRA3052

Cell Line Designation: A549 GAC23 C127

Date Sample Received: July 5, 2016  
Report Date: July 7, 2016

**Methodology:** Seventeen short tandem repeat (STR) loci plus the gender determining locus, Amelogenin, were amplified using the commercially available PowerPlex® 18D Kit from Promega. The cell line sample was processed using the ABI Prism® 3500Xi Genetic Analyzer. Data were analyzed using GeneMapper® ID-X v1.2 software (Applied Biosystems). Appropriate positive and negative controls were run and confirmed for each sample submitted.

**Data Interpretation** Cell lines were authenticated using Short Tandem Repeat (STR) analysis as described in 2012 in ANSI Standard (ASN-0002) *Authentication of Human Cell Lines: Standardization of STR Profiling* by the ATCC Standards Development Organization (SDO) and in Capes-Davis et al., *Match criteria for human cell line authentication: Where do we draw the line?* Int. J. Cancer. 2012 Nov 8. doi: 10.1002/ijc.27931

ATCC performs STR Profiling following ISO 9001:2008 and ISO/IEC 17025:2005 quality standards. There are no warranties with respect to the services or results supplied, express or implied, including, without limitation, any implied warranty of merchantability or fitness for a particular purpose. Neither ATCC nor Promega is liable for any damages or injuries resulting from receipt and/or improper, inappropriate, negligent or other wrongful use of the test results supplied, and/or from misidentification, misrepresentation, or lack of accuracy of those results. Your exclusive remedy against ATCC, Promega and those supplying materials used in the services for any losses or damage of any kind whatsoever, whether in contract, tort, or otherwise, shall be, at Promega's option, refund of the fee paid for such service or repeat of the service.

The ATCC trademark and trade name, any and all ATCC catalog numbers are trademarks of the American Type Culture Collection. PowerPlex is a registered trademark of Promega Corporation. Applied Biosystems, ABI Prism and GeneMapper are registered trademarks of Life Technologies Corporation.

Technical questions? Ordering questions?  
 ATCC Technical Support 800-438-6597 or 703-385-2700  
 (800) 638-6597 / +1 703-385-2700 Fax 703-385-2750  
 STRTechSupport@atcc.org Email: STRtesting@atcc.org

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Number of shared alleles between query sample and database profile: 15  
 Total number of alleles in the database profile: 15  
 Percent match between the submitted sample and the database profile: 100  
*The allele match algorithm compares the 8 core loci plus amelogenin only, even though alleles from all loci will be reported when available.*

**NOTE:** Loci highlighted in grey (8 core STR loci plus Amelogenin) can be made public to verify cell identity. In order to protect the identity of the donor, please do not publish the allele calls from all the STR loci tested. Electropherograms showing raw data are attached.

**Explanation of Test Results**  
 Cell lines with ≥ 80% match are considered to be related; i.e., derived from a common ancestry. Cell lines with between a 55% to 80% match require further profiling for authentication of relatedness.

The submitted sample profile is human, but not a match for any profile in the ATCC STR database  
 The submitted profile is an exact match for the following ATCC human cell line(s) in the ATCC STR database (8 core loci plus Amelogenin): CCL-185 (A549)  
 The submitted profile is similar to the following ATCC human cell line(s):

**Additional Comments:**  
 Submitted sample, STRA3052 (A549 GAC23 C127), is an exact match to ATCC cell line CCL-185 (A549).

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e-Signature, Technician: Isander 07/07/16

e-Signature, Reviewer: kkindig 07/07/16

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**Cell Line Authentication Service**  
STR Profile Report

FTA Barcode: STRA3052  
ATCC Sales Order: SOJ25378

Addendum: Comparative Output from the ATCC STR Profile Database

% Match	ATCC® Cell. No.	Designation	D5S818	D13S317	D7S820	D16S539	vWA	TH01	AMEL	TPOX	CSF1PO
100	CC-165	A549	11	11	8,11	11,12	14	8,9,3	X,Y	8,11	10,12
100	STRA3052	A549 GAC23 C127	11	11	8,11	11,12	14	8,9,3	X,Y	8,11	10,12

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**Definitions of terms used in this report:**

**Peak Area Difference (PAD):**  
Refers to a heterozygous peak imbalance. Two alleles at a single locus should amplify in a similar manner, and therefore produce peaks of similar height and area. Peaks which are above threshold (50 fls) but are not of similar area, within 50% of each other, are referred to as a PAD. Due to their nature cell lines do not amplify in the same manner as a sample taken from a fresh buccal swab. PAD is far more common in cell line samples.

**Stutter:**  
A stutter peak is a small peak which occurs immediately before the true peak. It is defined as being a single repeat unit smaller than the true peak. The stutter peak should be less than 15% of the true peak. The stutter is caused by the polymerase.

**+4 Peak:**  
A +4 is similar to a stutter but occurs immediately after the true peak. A stutter peak should be less than 5% for a homozygous and 10% for a heterozygous.

**Below Threshold Peak(s):**  
Cell lines can produce unusual profiles and occasionally a peak will amplify poorly and be below threshold. Where we find a below threshold peak which we believe is valid we indicate it as a below threshold peak. Our cell line analysis criteria, Homozygous and Heterozygous peaks must be equal to or above the set height threshold for it to be considered a true peak.

**Ladder/ Off Ladder Peak(s):**  
The allelic ladder consists of most or all known alleles in the population and allows for precise assignment of alleles. Those which do not align are termed 'off ladder'.

**Artifact:**  
A non-allelic product of the amplification process, an anomaly of the detection process, or a by-product of primer synthesis

**Pull-up:**  
A term used to describe when signal from one dye color channel produces artificial peaks in another, usually adjacent, color.

**Spike:**  
An extraneous peak resulting from dust, dried polymer, an air bubble, or an electrical surge.

**Dye blob:**  
Free dye not coupled to primer that can be injected into the capillary. (A known and documented dye blob is often found at the D3S1358 locus.)

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